



Desensitize with GLUMA®




Bond with iBOND®



Fill with Venus®



Cure with Translux®



**GLUMA® Desensitizer & GLUMA® Desensitizer PowerGel**  
Under every restoration.

To learn more, download  
the complimentary eBook:

**GlumaEbook.com**

Giving a hand to oral health.



**KULZER**  
MITSUI CHEMICALS GROUP

# GLUMA® Desensitizer

## GLUMA® Desensitizer PowerGel

Two solutions to cover all of your sensitivity needs



### GLUMA Desensitizer

GLUMA Desensitizer has been in the market for more than 20 years and has been used in more than 50 million restorations worldwide to stop and prevent hypersensitivity. One drop of Gluma is all you need and the result is fast and effective. No mixing, no curing and no repetitive steps.




### GLUMA Desensitizer PowerGel

With the same proven efficacy as GLUMA Desensitizer, GLUMA Desensitizer PowerGel is a unique one-step gel formula desensitizer that allows for greater control and accuracy during application. In addition, its unique green color indicator provides visual ease of placement and rinse-clean capability.

### How it works:

- **Desensitize:** GLUMA is the only desensitizer proven to penetrate exposed dentinal tubules up to 200µ. The result is the formation of multiple layers of protein septas that prevent inter-tubular fluid movements due to osmotic changes.
- **Inhibit bacteria:** GLUMA provides a hermetic seal that acts as a microbial barrier, inhibiting bacterial growth.
- **Re-wet:** GLUMA resurrects collapsed collagenous fibers, improving the bond strength of many adhesives.



Dentinal tubule before application of GLUMA



Protein septas following application of GLUMA

GLUMA® Desensitizer	Order Number
GLUMA Desensitizer 5ml	65872354
GLUMA Desensitizer Clinic Pack 3 x 5ml	66018221
GLUMA Desensitizer Single Dose 40 x 0.075ml	66001854
Application Brushes 2 Pack; 50 Pkg	66064661

GLUMA® Desensitizer PowerGel	Order Number
GLUMA PowerGel Syringe Kit 4 x 1gm	66043451
GLUMA PowerGel Brush Cannulas 60	66044334

**FIRST  
DESENSITIZER  
ON THE MARKET  
1987**



**DENTAL ADVISOR**  
1997 • 2000  
2008 (Top Desensitizer)



**REALITY'S CHOICE  
FIVE STAR AWARD**  
2001 • 2002 • 2003  
2005 • 2006 • 2007  
2009 • 2011



**TOWNIE CHOICE  
AWARD**  
2003 • 2004 • 2005  
2006 • 2007 • 2008  
2009 • 2010 • 2011

**CLINICIAN'S  
REPORT**  
2009

Top Effective Product for Class II

**GLUMA POWER GEL**



**DENTAL ADVISOR**  
2012



**CPS FIVE  
DIAMOND AWARD**  
2011



**REALITY'S CHOICE  
FIVE STAR AWARD**  
2013