

# Artic

Strong. Beautiful. Economical.

**PALA®**

## Mould Comparison

### Upper anteriors

Classic™	Artic®
134	O 20
136	S 30
137	S 35
1H	T 40
2D	O 10
3D	TO 05
2N	T 15
A24	TO 10
A25	T 20
A26	T 35
3M	S 05
3N	TO 20
3P	O 30
263	R 20
264	T 30
266	R 30
267	R 40

### Lower anteriors

Classic™	Artic®
2C	L 10
3M	L 05
2D	L 20
2N	L 15
3D	L 40
3N	L 25
2E	L 30
46	L 45
4H	L 55
3P	L 50

### Posteriors

Angulation	Classic™	Artic®
0°	29M	0°-30 U/L
0°	31M	0°-32 U/L
10°	F30	10°-30 U/L
10°	F32	10°-32 U/L
20°	29M	20°-30 U/L
20°	31M	20°-32 U/L
20°	33M	20°-34 U/L
20°	34M	20°-36 U/L
30°	30M	30°-30 U/L
30°	32M	30°-32 U/L
30°	34M	30°-34 U/L



**KULZER**  
MITSUI CHEMICALS GROUP

Classic™ is a trademark of Dentsply International Inc. Artic® is a trademark of Kulzer, LLC.