

NEW



## Introducing Delara®

The versatile mid-class tooth line with lifelike aesthetics.

Giving a hand to oral health.



**KULZER**  
MITSUI CHEMICALS GROUP



## Form, function and aesthetics.

Designed **by dental professionals, for dental professionals** using INCOMP and CAD/CAM technology, Delara offers these advantages:

### Beautiful, Lifelike Aesthetics

- Modern, bright shading (3 layer anteriors)
- Dynamic light refraction
- Vivid, enhanced surface structure
- Anatomical forms

### Easy Handling for Enhanced Productivity

- Polishing buffer
- Basal grinding reserve
- Optimized incisal edges
- Retains shape and layering even after grinding

### Versatility

- Optimized portfolio of shades and moulds with full range of options
- Various setup possibilities (tooth-to-tooth to cross bite)
- CAD/CAM production for easy integration into digital workflows

### Value and Performance

- Excellent functionality with Kulzer quality and durability
- Great price to performance ratio for mid-priced segment

*“Delara anterior teeth look so natural and beautiful in the mouth thanks to the unique surface structure. In addition, they are easy to set up, especially the posterior teeth. With Delara you have the freedom of set-up with almost no grinding necessary.”*



**Ilkka Garaisi**

As denturist (clinical dental technician) he has over 25 years of experience and treated over 7.000 patients with full dentures. He is president of EHT-federation in Finland, teacher of denturists, and representative of Finland in the IFD.



# Delara Anteriors

## SMALL

<b>T 42</b> ▽		42,4 9,6 8,0	<b>S 44</b> □		44,5 10,4 8,3
<b>O 43</b> ○		44,2 10,5 8,3	<b>T 44</b> ▽		45,1 10,5 8,2
<b>S 43</b> □		43,9 9,6 8,2	<b>R 45</b> □		45,0 10,6 8,2

## MEDIUM

<b>T 45</b> ▽		45,7 11,3 8,4	<b>O 47</b> ○		48,1 11,4 9,1
<b>T 46</b> ▽		46,5 10,7 8,7	<b>R 47</b> □		47,5 11,5 8,8
<b>R 46</b> □		46,5 10,1 8,8	<b>S 48</b> □		49,0 11,3 9,1

## LARGE

<b>R 48</b> □		48,8 13,0 9,3	<b>S 50</b> □		51,0 11,8 9,6
<b>T 49</b> ▽		49,6 12,2 9,1	<b>O 51</b> ○		52,2 12,4 9,8
<b>R 49</b> □		48,7 10,9 9,3	<b>R 52</b> □		52,1 13,6 9,8

**L 33**



33,6
9,4
5,2

**L 38**



36,6
10,5
5,8

**L 35**



34,6
9,5
5,4

**L 39**



38,3
10,8
5,9

**L 36**



35,6
9,1
5,5

## Delara Posteriors

**30**



29,5
9,1
31,4
8,1

**32**



31,0
9,5
33,1
8,8

**34**



33,0
10,2
35,3
9,3

**36**



34,8
10,9
37,3
9,9

### Mould numbers and possible combinations



Upper anteriors	
T 42	
T 44	
T 45	
T 46	
T 49	

Lower anteriors	
L 33	
L 33 + L 35	
L 33 + L 35 + L 36	
L 35 + L 36 + L 38	
L 38 + L 39	

Posteriors	
30 + 32	
30 + 32	
32 + 34	
32 + 34	
36	



O 43
O 47
O 51

L 33 + L 35 + L 36
L 38 + L 39
L 38 + L 39

30 + 32
34
36



R 45
R 46
R 47
R 48
R 49
R 52

L 33 + L 35
L 33 + L 35 + L 36
L 38 + L 39
L 36 + L 38
L 38 + L 39
L 38 + L 39

32
32 + 34
34
34
34 + 36
36




S 43
S 44
S 48
S 50


L 33 + L 35 + L 36
L 33 + L 35 + L 36 + L 38
L 38 + L 39
L 38 + L 39

30 + 32
30 + 32 + 34
34 + 36
36

### Examples



**T 42**  
42,4  
9,6  
8,0




**30 U**  
29,1

9,1



**L 33**  
33,6  
9,4  
5,2



**30 L**  
31,4

8,1

Also from Kulzer:

## Pala® Denture Teeth

### Perfection at your fingertips.

As one of the world's leading dental companies, Pala from Kulzer stands for over 80 years of prosthetic excellence. Pala is synonymous with high quality denture teeth, acrylics and first class accessories.

The development of these high quality, time-tested denture teeth is based on many years of experience, innovative technology and precision handcraft. Pala's large selection of mould designs provide the freedom and flexibility to select the perfect teeth for every patient's unique needs.

- **Mondial i** and **Mondial** denture teeth are made with the world's first and only dental acrylic—Flexecure®, a combination of unique microfillers and nanoparticles.
- **Artic's** multilayered structure allows exact matching of reflection curves in both natural and artificial teeth.
- **Artic Digital** offers CAD/CAM teeth at an economical price, specially designed to work with Kulzer's private processing economy digital dentures.

For more information, contact your  
Kulzer Sales Representative  
or visit [kulzerUS.com](http://kulzerUS.com).



**BRONZE**  
*Value*  
PALA Artic®



**SILVER**  
*Mid-class*  
Pala Delara®



**GOLD**  
*Premium*  
PALA Mondial®



**DIAMOND**  
*Ultra Premium*  
PALA Mondial® i



**Manufacturer**  
Kulzer GmbH  
Schlätterstraße 2  
88142 Wasserburg, Germany

Distributed by:  
**Kulzer, LLC**  
4315 S. Lafayette Blvd.  
South Bend, IN 46614  
(800) 431-1785